# Modell IR988

## Kontaktloses digitales Infrarot-Thermometer

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Guide and Manufacturer's Declaration -Electromagnetic Emissions		
Non-contact infrared thermometer is intended to be used in the following specified electromagnetic environments and the purchaser or user of the infrared thermometer shall ensure that it is used in this electromagnetic environment:		
Launch test	Compliance	Electromagnetic Environment-Guidance
Radio frequency emission GB 4824	1 Unit	Infrared thermometers use RF energy only for their internal functions. As a result, its RF emissions are low and the potential for interference with nearby electronic devices is low.
Radio frequency emission GB 4824	B type	Infrared thermometers are suitable for use in all installations, including domestic installations and direct connection to the public low-voltage supply network of domestic dwellings.
Harmonic emission GB 17625.1	Not Applicable	
Voltage fluctuation / flicker emission GB 17625.2	Not Applicable	

Recommended isolation distance between portable and mobile RF communications equipment and non-contact infrared thermometers			
Non-contact infrared thermometer is intended to be used in electromagnetic environments where RF radiation disturbances are controlled. According to the maximum rated output power of communication equipment, the purchaser or user may prevent electromagnetic interference by maintaining the minimum distance between the transmitter of portable and mobile radio frequency communication equipment and the non-contact infrared thermometer as recommended below.			
Rated maximum output power of transmitter/W	Isolation distance corresponding to different frequencies of the transmitter/m		
	150 kHz~80 MHz d=1.2√P	80 MHz~800 MHz d=1.2√P	800 MHz~2.5 GHz d=2.3√P
0.01	Not Applicable	0.12	0.23
0.1	Not Applicable	0.38	0.73
1	Not Applicable	1.2	2.3
10	Not Applicable	3.8	7.3
100	Not Applicable	12	23
For the transmitter rated maximum output power not listed in the table above, it is recommended that the isolation distance d, in meters (m), be determined by the formula in the corresponding transmitter frequency column, where P is the transmitter manufacturer's maximum output rated power, in watts (W). Note 1: At 80MHz and 800MHz, the formula of higher frequency range is used. Note 2: These guidelines may not be appropriate for all situations, where electromagnetic propagation is affected by the absorption and reflection of buildings, objects, and humans.			

Guidance and Manufacturer's Declaration-Electromagnetic Immunity			
Non-contact infrared thermometer is intended to be used in the following specified electromagnetic environments and the purchaser or user shall ensure that it is used in such electromagnetic environments:			
Immunity test	IEC60601 test level	Coincidence level	Electromagnetic Environment-Guidance
Radio frequency conduction GB/T 17626.6	3V (RMS) 150 kHz ~80 MHz	Not Applicable	Portable and Mobile RF Communication Devices shall not be used any part of the non-contact infrared thermometer including cables, closer than the recommended isolation distance. The distance shall be calculated by a formula corresponding to the transmitter frequency.  Recommended isolation distance: d=1.2√P 80 MHz~800 MHz d=2.3√P 800 MHz~2.5 GHz Wherein: P--Based on the transmitter manufacturer's maximum rated output power of the transmitter in Watt (W). d--Recommended isolation distance, in meter (m). The field strength of a fixed RF transmitter is determined by surveying the electromagnetic field, which shall be lower than the coincidence level in each frequency range. Interference may occur near devices marked with the following symbols:
Radio frequency radiation GB/T 17626.3 GB/T 17626.8	3V/m 80 MHz ~2.5 GHz	3V/m	
Note 1: At 80MHz and 800MHz, the formula for the higher frequency band is used. Note 2: These guidelines may not be appropriate for all situations, where electromagnetic propagation is affected by the absorption and reflection of buildings, objects, and humans.  a. The field strength of fixed transmitters, such as the field strength of wireless (cellular / cordless) telephones and ground mobile radio base stations, amateur radios, AM and FM radio broadcasts, and television broadcasts cannot be accurately predicted theoretically. In order to evaluate the electromagnetic environment of the fixed RF transmitter, the survey of electromagnetic field shall be considered. If the measured field strength of the non-contact infrared thermometer is higher than the applicable RF compliance level, the non-contact infrared thermometer shall be observed to verify its normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or positioning the non-contact infrared thermometer. b. The field strength shall be less than 3 V/m throughout the 150kHz ~ 80MHz frequency range.			

Guidance and Manufacturer's Declaration-Electromagnetic Immunity			
Non-contact infrared thermometer is intended to be used in the following specified electromagnetic environments and the purchaser or user shall ensure that it is used in such electromagnetic environments:			
Immunity test	IEC60601 test level	Coincidence level	Electromagnetic Environment-Guidance
Electrostatic discharge GB/T 17626.2	±6kV contact discharge ±8kV air discharge	±6kV contact discharge ±8kV air discharge	The floor shall be wood, concrete or tile, and if the floor is covered with synthetic material, the relative humidity shall be at least 30%.
Electrical fast transient burst GB/T 17626.4	±2kV to power line ±1kV to input / output line	Not Applicable	Not Applicable
Surge GB/T 17626.5	±1kV differential mode voltage ±2kV common mode voltage	Not Applicable	Not Applicable
Voltage sags, short interruptions and voltage changes on power input lines GB/T 17626.11	<5%U <sub>n</sub> , lasting for 0.5 cycles (On U <sub>n</sub> , >95% sag) 40% U <sub>n</sub> , lasting for 5 cycles (On U <sub>n</sub> , 60% sag) 70% U <sub>n</sub> , lasting for 25 cycles (On U <sub>n</sub> , 30% sag) <5% U <sub>n</sub> , lasting for 5s (On U <sub>n</sub> , >95% sag)	Not Applicable	Not Applicable
Power frequency magnetic field (150kHz) GB/T 17626.8	3A/m	3A/m,50/60Hz	The power frequency magnetic field shall have the power frequency magnetic field level characteristic of the typical place in the typical commercial or hospital environment.
Note: Ur refers to the AC network voltage before the test voltage is applied.			

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